

# Order Form

## Twin Block-Tool

An innovation from specialists for specialists



- A tool to simplify the manufacture of the Twin Block
- Exact 70 ° or 45 ° angel
- Parallel angels on the left and right block
- Saves time when spraying and trimming
- Individually adjusted, long lasting, economical, easy-care

I order a Twin Block-Tool for 99,- GBP or 129,95 EUR or 149,95 USD

**Instructions:**

70 ° by class II

45 ° by extreme class II

- Bend all wire on lower
- Lay the twin block-tool so that the magnets are on the distal third of the decidua or permanent five
- Check the upper plane of the twin block is parallel to the occlusal plane
- Press magnets gently onto the teeth thus close any gaps
- Wax with sticky wax the tool safely on the model
- Build wax wall mesial to the tool making the spraying of the bite planes easier.
- During the trimming process try and leave the inclined planes in their original condition
- Correct the height of the lower block by gently closing the Fixator/articulator while lowering the hite
- Trim and pumice the lower appliance
- Use ISOFOLAN from Scheu Dental to separate the upper and lower appliance while spraying the upper

Congratulations you have just made an accurate and correct Twin Block!

Hotline: 0049 6081 942131 or [team@reakfo.com](mailto:team@reakfo.com)

Contact ReakFO for technician courses on the Twin Block

Address, Stamp, Date, Signature

ReakFO Fachlaboratorium für Kieferorthopädie GmbH  
 In der Mark 53, 61273 Wehrheim, D - Deutschland  
 Tel: 0049 6081-942131  
**Fax: 0049 6081-942132**  
 Email: [team@reakfo.com](mailto:team@reakfo.com)

